## User's Manual of Face Shield


Face Shield is non-sterile disposable anti-fog Face Shield.

1. Brand: Refine
2. Product name: Face Shield
3. Material: Model A: PET 30% (Shield) PC 70% (Headwear),  
Model B: PET 80% PU sponge 20%.
4. Model: A, B (see packaging label for the detail model.)
5. Composition: composed of protective cover and fixing device made of polymer material.  
Non-sterile, single use.
6. Intended use : Face Shield is common eye protection equipment, which are used for protecting against droplets and splashes of liquids. Recommended industry: Facility Sanitation, Food Processing, Food Safety.
7. Product instruction: The product and the protective packaging should be checked before using. Stop using it if there is any damage.
8. Precautions: a. One-time use only. b. Use in accordance with the instructions. c. Never use it when damages are found. d. Stay away from chemicals.  
e. The eye protector is not intended to protect against high speed particles.
9. Storage: Please store in clean, dry and ventilated indoor place, with relative humidity (10%-93%), Atmospheric pressure for storage (70kPa-106kPa), temperature (-20°C- +40°C), avoid corrosive gas.

- Warning:** a. Materials which may come into contact with the wearer's skin could cause allergic reactions to susceptible individuals; b. Scratched or damaged oculars should be replaced;  
c. Eye-protectors against high speed particles worn over standard ophthalmic spectacles may transmit impacts, thus creating a hazard to the wearer.  
d. These protectors are intended for indoor use where no optical radiation hazards exist.


## Installation instructions (Model A)

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2. Please remove the protective membrane on both sides of the shield, then fix the shield on the frame by buckling up 4 buttons IN ORDER.

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## EU DECLARATION OF CONFORMITY

1. PPE (product, type, batch or serial number):

**Face Shield, A, Category II, AZ6A1**

2. Name and address of the manufacturer and, where applicable, his authorized representative:

**Manufacturer: Guilin Refine Medical Instrument Co.,LTD.**

**Address of the manufacturer: No.8-3, Information Industrial Park, High-Tech Zone, Qixing District, Guilin, Guangxi, 541004, P.R.China**

**Authorised representative: MedPath GmbH**

**Address of the Authorised representative: Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany**

3. This declaration of conformity is issued under the sole responsibility of the manufacturer:

**Manufacturer: Guilin Refine Medical Instrument Co.,LTD.**

**No.8-3, Information Industrial Park, High-Tech Zone, Qixing District, Guilin, Guangxi, 541004, P.R.China**

4. Object of the declaration (Identification of PPE allowing traceability; where necessary for the identification of the PPE a colour image of sufficient clarity may be included):

A:

Traceability Labeling:

**Face Shield**  
Brand: Refine  
Model: A  
Lot No.: **AZ6A1**  
Production Date: June 2020

5. The object of the declaration described in point 4 is in conformity with the relevant Union harmonization legislation: Regulation (EU) 2016/425

6. References to the relevant harmonized standards used, including the date of the standard, or references to the other technical specifications; including the date of the specification, in relation to which conformity is declared:

Harmonised Performance Standard No(s): **EN 166:2001**

Technical specification No(s): **RF-MIA-T001**

Test Reports: **C80202008R002**

7. Where applicable, the notified body UL International (Netherlands) B.V. (European Notified Body No. 2821) performed the EU type-examination Module B and issued the EU type-examination certificate : 2821-PPE-0002.

Signed for and on behalf of Guilin Refine Medical Instrument Co.,LTD.

(place and date of issue): Guilin , 2020-06-07

Name: Jordan Chen, Title: Management representative (signature):

## EU DECLARATION OF CONFORMITY

1. PPE (product, type, batch or serial number):

**Face Shield, B, Category II, AZ6B1**

2. Name and address of the manufacturer and, where applicable, his authorized representative:

**Manufacturer: Guilin Refine Medical Instrument Co.,LTD.**

**Address of the manufacturer: No.8-3, Information Industrial Park, High-Tech Zone, Qixing District, Guilin, Guangxi, 541004, P.R.China**

**Authorised representative: MedPath GmbH**

**Address of the Authorised representative: Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany**

3. This declaration of conformity is issued under the sole responsibility of the manufacturer:

**Manufacturer: Guilin Refine Medical Instrument Co.,LTD.**

**No.8-3, Information Industrial Park, High-Tech Zone, Qixing District, Guilin, Guangxi, 541004, P.R.China**

4. Object of the declaration (Identification of PPE allowing traceability; where necessary for the identification of the PPE a colour image of sufficient clarity may be included):

B:

Traceability Labeling:

**Face Shield**

Brand: Refine

Model: B

Lot No.: **AZ6B1**

Production Date: June 2020

5. The object of the declaration described in point 4 is in conformity with the relevant Union harmonization legislation: Regulation (EU) 2016/425

6. References to the relevant harmonized standards used, including the date of the standard, or references to the other technical specifications; including the date of the specification, in relation to which conformity is declared:

Harmonised Performance Standard No(s): **EN 166:2001**

Technical specification No(s): **RF-MIB-T001**

Test Reports: **C80282070**

7. Where applicable, the notified body UL International (Netherlands) B.V. (European Notified Body No. 2821) performed the EU type-examination Module B and issued the EU type-examination certificate : 2821-PPE-0003.

Signed for and on behalf of Guilin Refine Medical Instrument Co.,LTD.

(place and date of issue): **Guilin , 2020-06-07**

Name: Jordan Chen, Title: Management representative (signature):

*Jordan Chen*